

Enterotube[®] differentiates

Perform 9 Tests on 8 Separate Standard Media

and identify 17 members of the Enterobacteriaceae family in one tube. The single enclosed system with its own inoculating needle avoids error and contamination. Incubate only 18 to 24 hours for most tests and easily interpret the diagnostic color changes in each medium.

Shorten Diagnosis 24 Hours and Economize

Inoculate directly from primary culture and save a day. Obtain more information for differentiating and identifying Enterobacteriaceae than possible with TSI within the same time period. Small colonies are sufficient to inoculate ENTEROTUBE plus additional tests. Eliminate media preparation, glassware washing and sterilization, and save time and space.

